

PRODUCT BROCHURE

KERN 317-08-2022i

Date: 15.03.2025

ENGLISH



Individual weight, OIML Class E2, 200 g, Knob, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	317-08-2022i
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Trade parts list element
EAN code (Model)	4045761245341
Predecessor Model	317-08; 317-08-2021a
Successor Model	317-08
Customs tariff number	84239090

Technical data - Primary

Nominal value	200 g
Tolerance (OIML)	0,3 mg
Case aluminium	317-080-600
Case plastic	317-080-400
Case wood	317-080-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

OIML class	E2
OIML conform	yes
Density (kg/m ³)	7950
Volume (cm ³)	25

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	28×28×47 mm
Gross weight	0.20 kg
Net weight	0.20 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-338
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
317-080-100	Wooden Box for Set of Weights KERN 317-080-100
317-080-400	Plastic Weight Box KERN 317-080-400
317-080-600	Aluminium Protective Box KERN 317-080-600
315-242	Tweezers KERN 315-242

SERVICES

Model	Description
952-338	Verification KERN 952-338
962-338	DAkkS Calibration KERN 962-338
962-338R	DAkkS Recalibration KERN 962-338R
920-33	Marking by etching KERN 920-33
921-338	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-338
920-331	Marking by etching (last digits of the serial number) KERN 920-331