

PRODUCT BROCHURE

KERN 313-07

Date: 09.03.2026

ENGLISH



Weight set, E2, 1 mg - 2 kg, Knob, stainless steel polished, in wooden case


CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information



Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	305×190×140 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-306
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-306</u>	DAkkS Calibration KERN 962-306
<u>962-306R</u>	DAkkS Recalibration KERN 962-306R
<u>921-306</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-306
<u>920-707</u>	Marking (Last digits of serial number) KERN 920-707

SPARE PARTS

Model	Description
<u>313-070-100</u>	Box of Weights KERN 313-070-100
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270