


PRODUCT BROCHURE

KERN 346-065-009

Date: 17.04.2025

ENGLISH



Weight container for rectangular weights, calibrated to OIML class M1

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- Preconfigured weight container for testing high-load floor scales, pallet scales, pallet truck scales, crane scales, etc. This can also be used for storing the weights.
- This means the weight container and the weights can be placed on the balance in one go, saving you time and money
- The weight container is adjusted to OIML accuracy class M1. Other OIML accuracy classes are also available, please ask
- Weight of the weight container, OIML class M1: 60 kg
- Possible equipment, rectangular weights, OIML class M1: 8×50 kg + 2×20 kg
- Maximum total weight (weight container incl. weights): 500 kg
- Tolerance ± 3 g

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	346-065-009
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1225 107
Customs tariff number	84239090

Technical data - Primary

Nominal value	60 kg
Tolerance (OIML)	3 g

Technical data - Construction

Dimensions	1000×500×500 mm
------------	-----------------

Material stainless steel

Technical data - Measuring system

Nominal value	60 kg
OIML class	M1
OIML conform	yes
Tolerance (OIML)	3 mg

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	1000×500×500 mm
Gross weight	60 kg
Net weight	60 kg
Shipping method	Freight
Packaging component - by weight - plastics	0 g
Delivery time	2 d

Optional Services

DAkkS Calibration 962-646

FUNCTIONS

Standard



SERVICES

Model	Description
962-646	DAkkS Calibration KERN 962-646
962-646R	DAkkS Recalibration KERN 962-646R

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium
346-06	5 kg	250 mg	346-060-600
346-07	10 kg	500 mg	346-070-600
346-08	20 kg	1 g	346-080-600
346-09	50 kg	2,5 g	346-090-600
346-022-005	20 kg	1 g	-
346-042-008	40 kg	1,5 g	-
346-052-010	50 kg	2,5 g	-
346-055-004	50 kg	2,5 g	-
346-055-009	50 kg	2,5 g	-
346-065-009	60 kg	3 g	-
346-86	5 kg	250 mg	346-060-600
346-87	10 kg	500 mg	346-070-600
346-77	10 kg	500 mg	346-070-600
346-76	5 kg	250 mg	346-060-600
346-79-2023a	50 kg	2,5 g	-
346-88	20 kg	1 g	346-080-600
346-78	20 kg	1 g	346-080-600
346-89	50 kg	2,5 g	-

Model name	Nominal value	Tolerance (OIML)	Case aluminium
346-81	100 kg	5 g	-
346-82	200 kg	10 g	-
346-83	500 kg	25 g	-
346-84	1000 kg	50 g	-
346-79	50 kg	2,5 g	346-090-600