

PRODUCT BROCHURE

KERN 347-117-100

Date: 11.03.2026

ENGLISH



Supporting bar, OIML Class M1, 1 kg, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- for fixing slotted weights
- A precise curvature at the lower end of the beam bar ensures that the slotted weights rest securely. This design prevents the weights from slipping and ensures stability – even in slight inclined positions
- Own weight beam bar 1 kg

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-117-100
Model-Serie	347-7
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1366756
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 kg
Tolerance (OIML)	50 mg

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

OIML class	M1
OIML conform	yes
Density (kg/m ³)	8000
Volume (cm ³)	119,0476

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	55×55×384 mm
Gross weight	1.00 kg
Net weight	1.00 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-641
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>347-037</u>	Test Weight KERN 347-037
<u>347-047</u>	Test Weight KERN 347-047
<u>347-057</u>	Test Weight KERN 347-057
<u>347-067</u>	Test Weight KERN 347-067
<u>347-077</u>	Test Weight KERN 347-077
<u>347-087</u>	Test Weight KERN 347-087
<u>347-097</u>	Test Weight KERN 347-097
<u>347-117</u>	Test Weight KERN 347-117
<u>347-127</u>	Test Weight KERN 347-127
<u>347-137</u>	Test Weight KERN 347-137
<u>347-147</u>	Test Weight KERN 347-147

SERVICES

Model	Description
<u>962-641</u>	DAkkS Calibration KERN 962-641
<u>962-641R</u>	DAkkS Recalibration KERN 962-641R