

PRODUCT BROCHURE

KERN 308-426

Date: 10.03.2026

ENGLISH



Set of milligram weights, E1, 1 mg - 500 mg, Wire, stainless steel polished, in aluminium case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	308-426
Model-Serie	308
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main article
EAN code (Model)	4045761223608
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 500 mg
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 mg - 500 mg
OIML class	E1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	165×158×37 mm
Gross weight	0.35 kg
Net weight	0.35 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-250
-------------------	---------

FUNCTIONS

Standard



Option



ENGLISH

SERVICES

Model	Description
<u>962-250</u>	DAkkS Calibration KERN 962-250
<u>962-250R</u>	DAkkS Recalibration KERN 962-250R
<u>920-701</u>	Marking (Last digits of serial number) KERN 920-701

SPARE PARTS

Model	Description
<u>313-010-600</u>	Aluminium Protective Box KERN 313-010-600
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270