

PRODUCT BROCHURE

KERN 366-96

Date: 11.03.2026

ENGLISH



Individual weight, OIML Class M3, 5 kg, Knob, cast iron lacquered

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	366-96
Model-Serie	366
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1482 104
Predecessor Model	366-96-2024a
Customs tariff number	84239090
Country of Origin	IN

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	2,5 g

Technical data - Construction

Material	cast iron lacquered
----------	---------------------

Technical data - Measuring system

Nominal value	5 kg
OIML class	M3
OIML conform	yes
Tolerance (OIML)	2,5 g
Density (kg/m ³)	7100

Volume (cm³) 704,2254

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	94×94×146 mm
Gross weight	5 kg
Net weight	5 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	no
-----------------------	----

Optional Services

DAkkS Calibration	962-643
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
317-280	Gloves KERN 317-280

SERVICES

Model	Description
962-643	DAkkS Calibration KERN 962-643
962-643R	DAkkS Recalibration KERN 962-643R