

KERN 367-01

Individual weight, OIML Class M3, 1 g, Knob, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	367-01
Model-Serie	367-01
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761257207
Predecessor Model	367-41
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g
Tolerance (OIML)	10 mg
Case aluminium	317-010-600
Case plastic	347-030-400
Case wood	337-010-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	1 g
OIML class	M3
OIML conform	Yes

Functions

Standard

Tolerance (OIML) 10 mg

Density (kg/m³) 7.950

Volume (cm³) 0,119

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 6×6×6 mm

Gross weight 0,05 kg

Net weight 0,05 kg

Shipping method Parcel service

Packaging component - by weight - plastics 0 g

Delivery time 1 d

Technical data - Verification

Verification possible No

Optional Services

DAkkS Calibration 962-631

KERN 367-01

Individual weight, OIML Class M3, 1 g, Knob, stainless steel fine turned



Option

