

KERN 308-32

Milligram weights, E1, 2 mg, Wire, stainless steel polished



Article codes, product groups & Status information

Model article-code	308-32
Model-Serie	308
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761480988
Predecessor Model	308-32-2.024a
Planned start date	2025-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	2 mg
Tolerance (OIML)	0,003 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	6×6×0,3 mm
Material	Stainless steel polished

Technical data - Measuring system

Nominal value	2 mg
OIML class	E1

Functions

Standard

OIML conform Yes

Tolerance (OIML) 0,003 mg

Density (kg/m³) 8.000

Volume (cm³) 0,0003

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 6×6×3 mm

Gross weight 0,05 kg

Net weight 0,05 kg

Shipping method Parcel service

Packaging component - by weight - plastics 0 g

Delivery time 1 d

Optional Services

DAkkS Calibration 962-252

KERN 308-32

Milligram weights, E1, 2 mg, Wire, stainless steel polished



Option

