

PRODUCT BROCHURE

KERN 317-11

Date: 09.03.2026

ENGLISH



Individual weight, OIML Class E2, 1 kg, Knob, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class E2, 1 kg, Knob, stainless steel polished (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section
- For the verification you need a box. Please add a box for single weights.

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	317-11
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1003705
Successor Model	317-11-2022i
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 kg
Tolerance (OIML)	1,6 mg
Case aluminium	317-110-600
Case plastic	317-110-400
Case wood	317-110-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	1 kg
OIML class	E2
OIML conform	yes
Tolerance (OIML)	1,6 mg
Density (kg/m ³)	7950
Volume (cm ³)	125

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	48×48×82 mm
Gross weight	1.00 kg
Net weight	1.00 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible yes

Optional Services

DAkkS Calibration 962-341

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>317-280</u>	Gloves KERN 317-280
<u>317-110-100</u>	Wooden Box for Set of Weights KERN 317-110-100
<u>317-110-400</u>	Plastic Weight Box KERN 317-110-400
<u>317-110-600</u>	Aluminium Protective Box KERN 317-110-600

SERVICES

Model	Description
<u>962-341</u>	DAkkS Calibration KERN 962-341
<u>962-341R</u>	DAkkS Recalibration KERN 962-341R
<u>921-341</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-341