

## PRODUCT BROCHURE

---

# OPTICS OZB-A5203

Date: 30.04.2025

ENGLISH



Individuality, variety and flexible working through our modular construction system ? PREMIUM-Universal stands

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

Stereomicroscope Stand OZB-A5203

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                 |
|-----------------------|---------------------------------|
| Model article-code    | OZB-A5203                       |
| Model-Serie           | OZB-UP                          |
| Assortment            | Optics                          |
| Product Type          | Microscopes                     |
| Product group         | Stereo microscopes              |
| Product usage type    | Spare part;Accessories/Services |
| EAN code (Model)      | 4045761169593                   |
| Customs tariff number | 90119000                        |

### Technical data - Construction

|                    |                                                              |
|--------------------|--------------------------------------------------------------|
| Stand type         | universal (ball-bear<br>double telescopic arm<br>with plate) |
| Stand area         | 254×254 mm                                                   |
| Holder diameter    | 76 mm                                                        |
| Dimensions (W×D×H) | 254×545×515 mm                                               |

### Technical data - Environmental conditions

|                  |       |
|------------------|-------|
| Storage Temp Min | -5 °C |
| Stotage Temp Max | 40 °C |

### Technical data - Packing & Shipping

|                                             |                |
|---------------------------------------------|----------------|
| Dimensions packaging (W×D×H)                | 680×330×260 mm |
| Gross weight                                | 24 kg          |
| Net weight                                  | 22 kg          |
| Shipping method                             | Parcel service |
| Packaging component - by weight - paper     | 0 g            |
| Packaging component - by weight - cardboard | 699.761 g      |
| Packaging component - by weight - plastics  | 23.8104 g      |
| Packaging component - by weight - styrofoam | 599.228 g      |
| Packaging component - by weight - foam      | 0 g            |
| Delivery time                               | 1 d            |

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|



## ACCESSORIES

| Model     | Description                                               |
|-----------|-----------------------------------------------------------|
| OZB-A5301 | Stereo Microscope Modular System – holders KERN OZB-A5301 |
| OZB-A5306 | Stereo Microscope Modular System – holders KERN OZB-A5306 |

## SPARE PARTS

| Model     | Description                   |
|-----------|-------------------------------|
| OZB-A5806 | Replacement lever OZB-A5806   |
| OZB-A5290 | Plate + Pillar KERN OZB-A5290 |
| OZB-A5291 | Plate KERN OZB-A5291          |