


KERN ODC 841

Specialists in microscopy for measurement, counting, documentation, archiving and image processing



Article codes, product groups & Status information

Model article-code	ODC 841
Model-Serie	ODC-84
Assortment	Optics
Product Type	Microscopes
Product group	Microscope cameras
Product usage type	Main article
EAN code (Model)	4045761235724
Planned start date	2020-1HJ
Customs tariff number	85258900

Technical data - Primary

Camera type	USB3 Microscope camera - C-Mount
Camera resolution	20 MP
Sensor type	Sony CMOS
Sensor size	1"

Technical data - Construction

Color capability	Colour
Camera mounting	C-Mount
Interfaces	USB 3.0
Operating system	Win 10;Win 11;Win 7;Win 8
Applicable software	External (OXM 902)
Dimensions	68×68×46 mm

Technical data - Power supply

Power Plug	USB B
------------	-------

Protection Class

Class III - Low voltage

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	150×150×100 mm
Gross weight	0.65 kg
Net weight	0.50 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	91.5 g
Packaging component - by weight - plastics	500 mg
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - glass	0 g
Packaging component - by weight - wood	0 g
Packaging component - by weight - foam	72.5 g
Packaging component - by weight - metal	0 g
Delivery time	1 d

Technical data - Environmental conditions

KERN ODC 841

Specialists in microscopy for measurement, counting, documentation, archiving and image processing



Storage Temp Min -5 °C

Storage Temp Max 40 °C

Technical data - Verification

WEEE Yes

Functions

Standard

