

KERN 347-126

Individual weight, OIML Class M1, 2 kg, Hook weight, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-126
Model-Serie	347-126
Assortment	Scales
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761256255
Predecessor Model	347-526
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	2 kg
Tolerance (OIML)	100 mg
Case plastic	347-130-400

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	2 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	100 mg
Density (kg/m ³)	8.000

Functions

Standard

Volume (cm³) 238,0952

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	60×60×146 mm
Gross weight	2,0 kg
Net weight	2,0 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-642
-------------------	---------

KERN 347-126

Individual weight, OIML Class M1, 2 kg, Hook weight, stainless steel fine turned



Option

