

PRODUCT BROCHURE

KERN 357-13

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M2, 5 kg, Knob, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	357-13
Model-Serie	357-13
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761289048
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	800 mg
Case aluminium	317-130-600
Case plastic	347-130-400
Case wood	337-130-200

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	5 kg
OIML class	M2

OIML conform	yes
Tolerance (OIML)	800 mg
Density (kg/m ³)	7950
Volume (cm ³)	595,2381

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	100×100×165 mm
Gross weight	5 kg
Net weight	5 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-643
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>317-280</u>	Gloves KERN 317-280
<u>317-290</u>	Gloves KERN 317-290
<u>337-130-200</u>	Wooden Box for Set of Weights KERN 337-130-200
<u>347-130-400</u>	Plastic Weight Box KERN 347-130-400
<u>317-130-600</u>	Aluminium Protective Box KERN 317-130-600

SERVICES

Model	Description
<u>962-643</u>	DAkkS Calibration KERN 962-643
<u>962-643R</u>	DAkkS Recalibration KERN 962-643R