

PRODUCT BROCHURE

KERN 315-412-408

Date: 19.03.2025

ENGLISH



Safety Set 200 g (M1); 2000 g (M1)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



DESCRIPTION

- These “KERN Safety Sets” each contain one test weight for testing the sensitivity, i.e. the correct adjustment of your balance as well as a small test weight to test the lower end of the weighing range, the so-called minimum sample weight
- As an option, the “KERN Safety Set” has space for another test weight, for testing your balance at a weight which is relevant for you
- Useful accessories which have been selected to suit that particular “KERN Safety Set”, such as, for example, special gloves, tweezers, weight grips, brushes, etc. will assist you in handling your test weights properly
- Stored in the practical protective case next to your balance, you can check and ensure the high precision of your balance at any time

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	315-412-408
Model-Serie	315
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services

EAN code (Model) 4045761245587

Customs tariff number 84239090

Technical data - Primary

Nominal value 200 g & 2 kg

Case aluminium zet315-000-000

Weight set case in aluminium case

Technical data - Construction

Material stainless steel

Is weight set	yes
Number of weights in set	2
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	200 g & 2 kg
OIML class	M1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	255×135×140 mm
Gross weight	3 kg
Net weight	1.2 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	1.8 kg
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-600-600
-------------------	-------------

FUNCTIONS

Standard



1 DAY

SERVICES

Model	Description
962-600-600	DAkkS Calibration KERN 962-600-600

SPARE PARTS

Model	Description
315-242	Tweezers KERN 315-242

RELATED ARTICLES DATA

Model name	Nominal value	Case aluminium	Weight set case
315-406-403	5 g & 50 g	zet315-000-000	in aluminium case
315-412-407	100 g & 2 kg	zet315-000-000	in aluminium case
315-409-405	20 g & 500 g	zet315-000-000	in aluminium case
315-412-408	200 g & 2 kg	zet315-000-000	in aluminium case
315-407-404	10 g & 100 g	zet315-000-000	in aluminium case
315-408-404	10 g & 200 g	zet315-000-000	in aluminium case
315-408-405	20 g & 200 g	zet315-000-000	in aluminium case
315-409-406	50 g & 500 g	zet315-000-000	in aluminium case
315-411-407	100 g & 1 kg	zet315-000-000	in aluminium case
315-411-406	50 g & 1 kg	zet315-000-000	in aluminium case