

KERN MPB 300K100P

Personal floor scale with BMI function



Article codes, product groups & Status information

Model article-code	MPB 300K100P
Model-Serie	MPB
Assortment	KERN
Product Type	Medical products
Product group	Personal scales
Product usage type	Main article
EAN code (Model)	4045761097582
Planned start date	2008-1HJ
Customs tariff number	84231090
CE Label	Yes

Technical data - Primary

Weighing capacity	300 kg
Readability	0.1 kg

Technical data - Construction

Material weighing plate	Steel, lacquered
Dimension completely mounted (W×D×H)	315×450×1000 mm
Material casing	Steel, lacquered

Technical data - Display

Display height	25 mm
----------------	-------

Technical data - Power supply

Battery	6×1.5 V AA
Operating time	170 h

Mains adapter/adaptor included

CH;EURO

Input Voltage	220 V - 240 V AC, 50 Hz
---------------	-------------------------

Technical data - Measuring system

Calibration method	Adjusting with external weight
Linearity	± 0.2 kg
Readability	0.1 kg
Recommended adjusting weight	250 kg (M1)
Repeatability	0.1 kg
Weighing units	Kg
Weighing system	Strain gauge

Technical data - Environmental conditions

Maximum humidity	80 %
Maximum operating temperature	35 °C
Minimum ambient temperature	5 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	1175×410×330 mm
Gross weight	10 kg
Net weight	7 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	720 g
Packaging component - by weight - plastics	1.8 kg

KERN MPB 300K100P

Personal floor scale with BMI function



Packaging component - by weight - styrofoam 2880 g

Delivery time 1 d

Optional Services

DAkkS Calibration 963-129

Functions

Standard



Option

