

KERN TVL 500FHS71

The practical all-in-one package for rapid, simple testing



Article codes, product groups & Status information

Model article-code	TVL 500FHS71
Model-Serie	TVL
Assortment	Sauter
Product Type	Force measurement
Product group	Test stands (mechanic)
Product usage type	Main article
EAN code (Model)	4045761273641
Predecessor Model	AE 500-S02
Planned start date	2021-1HJ
Customs tariff number	90248000
CE Label	Yes

Technical data - Primary

Max. Force	500 N
Max. Travel distance	230 mm

Technical data - Measuring system

Max. Force	500 N
Readability	0.1 N
Tolerance (% from Max.)	0.5%
Maximum measured length	20 cm

Technical data - Environmental conditions

Maximum storage temperature	40 °C
Minimum storage temperature	-10 °C

Functions

Standard

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	570×275×300 mm
Gross weight	11 kg
Net weight	10 kg
Shipping method	Parcel service
UN number	NiMH (7 - not restricted as per A199)
Delivery time	1 d

KERN TVL 500FHS71

The practical all-in-one package for rapid, simple testing



Option

