

PRODUCT BROCHURE

KERN KUP-04

Date: 16.03.2025

ENGLISH



Interface adapter with cable. Interface with Ethernet
KUP, cable length 0,15 m

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- The KERN Universal Port (KUP) enables the exchange of data and control commands between balances and the network or PC, mobile end unit, printer, etc. using several interfaces: RS-232, USB device, Ethernet, Bluetooth and WiFi. In this way we ensure the perfect integration of your KERN balance into production or processing chains, which gives the benefit of a simple, consistent workflow
- In accordance with the „Plug and Play“ principle, the KUP for RS-232 or USB, Ethernet, Bluetooth or WiFi will only be inserted into the DB 15 socket which is integrated into the device. It is therefore possible to retrofit or replace interfaces conveniently without having to open the balance casing and no complicated installation. This option therefore offers a high level of flexibility at low additional cost

SPECIFICATIONS

Article codes, product groups & Status information

| | |
|-----------------------|------------------------------|
| Model article-code | KUP-04 |
| Model-Serie | KUP |
| Assortment | KERN |
| Product Type | General weighing accessories |
| Product group | Accessories (printer, etc.) |
| Product usage type | Accessories/DL |
| EAN code (Model) | 4045761253933 |
| Customs tariff number | 90160090 |

Packaging component - 32 g
by weight - styrofoam

Delivery time 1 d

FUNCTIONS

Standard



Technical data - Packing & Shipping

| | |
|---------------------------------------------|----------------|
| Dimensions packaging (W×D×H) | 225×145×40 mm |
| Gross weight | 0.15 kg |
| Net weight | 0.10 kg |
| Shipping method | Parcel service |
| Packaging component - by weight - cardboard | 32 g |