

## PRODUCT BROCHURE

---

# KERN 336-36

Date: 19.03.2025

ENGLISH



Individual weight, OIML Class F2, 5 kg, Block, stainless steel

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)

## DESCRIPTION

- Individual weight, OIML Class F2, 5 kg, Block, stainless steel
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section
- For the verification you need a box. Please add a box for single weights.

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 336-36                         |
| Model-Serie           | 336                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML F2 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761232808                  |
| Customs tariff number | 84239090                       |

### Technical data - Primary

|                  |               |
|------------------|---------------|
| Nominal value    | 5 kg          |
| Tolerance (OIML) | 80 mg         |
| Case aluminium   | 346-060-600   |
| Case plastic     | not available |
| Case wood        | not available |

### Technical data - Construction

|            |                 |
|------------|-----------------|
| Dimensions | 150×68×96 mm    |
| Material   | stainless steel |

### Technical data - Measuring system

|                  |       |
|------------------|-------|
| Nominal value    | 5 kg  |
| OIML class       | F2    |
| OIML conform     | yes   |
| Tolerance (OIML) | 80 mg |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 150×68×96 mm   |
| Gross weight                               | 5 kg           |
| Net weight                                 | 5 kg           |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

### Technical data - Verification

|                       |     |
|-----------------------|-----|
| Verification possible | yes |
|-----------------------|-----|

### Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-443 |
|-------------------|---------|

## FUNCTIONS

Standard



## ACCESSORIES

| Model       | Description                               |
|-------------|-------------------------------------------|
| 346-060-600 | Aluminium Protective Box KERN 346-060-600 |

## SERVICES

| Model    | Description                                                            |
|----------|------------------------------------------------------------------------|
| 962-443  | DAkkS Calibration KERN 962-443                                         |
| 962-443R | DAkkS Recalibration KERN 962-443R                                      |
| 920-30   | Marking with stroke number KERN 920-30                                 |
| 921-443  | Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-443                   |
| 920-301  | Marking with stroke number (Last digits of serial number) KERN 920-301 |
| 920-321  | Marking by etching (last digits of the serial number) KERN 920-321     |

## RELATED ARTICLES DATA

| Model name | Nominal value | Tolerance (OIML) | Case aluminium | Case plastic  | Case wood     |
|------------|---------------|------------------|----------------|---------------|---------------|
| 336-36     | 5 kg          | 80 mg            | 346-060-600    | not available | not available |
| 336-37     | 10 kg         | 160 mg           | 346-070-600    | not available | not available |
| 336-38     | 20 kg         | 300 mg           | 346-080-600    | not available | not available |
| 336-39     | 50 kg         | 800 mg           | 346-090-600    | not available | not available |