


Software EasyTouch

SET-46

Backup Scheduler



EasyTouch Backup Scheduler – data backup and restore

Features

- Note: Prerequisite for this set is the basic program SET-01 Base
- Backup Scheduler: Automate your data backups and reap the benefits of continuous, scheduled processes which keep your valuable information safe in the background
- Automatic data backups: Save time and money through automatic backups at defined intervals. This ensures consistent, reliable record keeping, reduces downtime and minimises the risk of data loss in the event of disruptions
- Backup file configuration: Users can define the name of the backup file individually, including a master name and a timestamp as a postfix. The backup directory and type of backup can also be individually adapted – whether it's the creation of a new file, or overwriting or expanding an existing file
- Backup interval: The interval for backups can be set in hours, several times or once a day, weekly, monthly and annually
- Backup messages: Every automatic data backup will be displayed in the backup directory with the file name as a link to the backup file
- Backup history: The progress and status of all backups are logged and can be traced at any time
- Restore: In the event of a fault, the last saved backup can be easily restored using the “Restore” function

STANDARD



OPTION

