

KERN 307-07

Individual weight, OIML Class E1, 100 g, Knob, stainless steel polished



Article codes, product groups & Status information

Model article-code	307-07
Model-Serie	307
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761314245
Predecessor Model	37-7-2.021a
Planned start date	2022-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	100 g
Tolerance (OIML)	0,05 mg
Case aluminium	317-070-600
Case plastic	317-070-400
Case wood	317-070-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	100 g
OIML class	E1
OIML conform	Yes

Functions

Standard

Tolerance (OIML)	0,05 mg
Density (kg/m ³)	8.000
Volume (cm ³)	12,5

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	22×22×38 mm
Gross weight	0,10 kg
Net weight	0,10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-237
-------------------	---------

KERN 307-07

Individual weight, OIML Class E1, 100 g, Knob, stainless steel polished



Option

