


PRODUCT BROCHURE

OPTICS OZB-A1213

Date: 30.04.2025

ENGLISH



Individuality, variety and flexible working through our modular construction system ? ECO-Universal stands

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Stereomicroscope Stand OZB-A1213

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OZB-A1213
Model-Serie	OZB-UE
Assortment	Optics
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761218314
Customs tariff number	90119000
CE Label	yes

Dimensions packaging (W×D×H)	490×420×290 mm
------------------------------	----------------

Gross weight	17 kg
--------------	-------

Net weight	15 kg
------------	-------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - cardboard	1 kg
---	------

Packaging component - by weight - plastics	1 kg
--	------

Delivery time	1 d
---------------	-----

FUNCTIONS

Standard



Technical data - Construction

Stand type	universal (ball-beared double telescopic arm with clamp)
Holder diameter	76 mm
Dimensions (W×D×H)	526×210×380 mm

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

ENGLISH

SPARE PARTS

Model	Description
OZB-A1902	Replacement lever OZB-A1902
OZB-A1903	Replacement lever OZB-A1903
OZB-A1906	Holder connecting rod OZB-A1906
OZB-A1901	Replacement lever OZB-A1901