

PRODUCT BROCHURE

KERN YKI-12

Date: 19.09.2025

ENGLISH



RS-232/WiFi-Converter

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- Converter for wireless connection of scales, force measuring instruments, etc. with RS-232 interface to WiFi-enabled devices such as tablets, laptops, smartphones or networks
- In this way, where there is a physical separation or a large distance between the devices, the weighing data can be easily transmitted wirelessly from a scale to one or more WiFi devices or to a network (corresponding software required on the WiFi device, not included with the delivery)
- Allows direct integration of measuring devices into most company networks, with automatic configuration using DHCP or statistical IP allocation. The current measurement data can be accessed from the entire network at any time
- The converter supports a wide range of configuration options via the web interface
- Compatible with KERN BalanceConnection software
- Included with delivery: RS-232/WLAN converter, mains adapter (EU)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YKI-12
Model-Serie	YKI
Assortment	KERN
Product Type	General weighing accessories

Product group	Accessories (printer, etc.)
Product usage type	Accessories/DL
EAN code (Model)	4045761483385
Predecessor Model	YKI-03
Customs tariff number	84239010
CE Label	yes

Technical data - Power supply

Mains adapter	Power supply unit
---------------	-------------------

Technical data - Packing & Shipping

Gross weight	0.40 kg
Net weight	0.20 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	194 g
Delivery time	1 d

FUNCTIONS

Standard



SPARE PARTS

Model	Description
YKA-11	Power supply KERN YKA-11