

KERN FOB 10K-3LM

Compact stainless steel bench scale with special checkweighing display for even more efficient working, verification optional



Article codes, product groups & Status information

Model article-code	FOB 10K-3LM
Model-Serie	FOB-LM
Assortment	Scales
Product Type	Food scales
Product group	Bench scales
Product usage type	Main article
EAN code (Model)	4045761214354
Successor Model	FFN 15K5IPM; FKS 15K-3M
Planned start date	2019-1HJ
Customs tariff number	84238129
Model Status	Yes
CE Label	Yes
Country of Origin	CN

Technical data - Primary

Weighing capacity	15.000 g
Readability	5 g
Verification value (e)	5 g
Minimum load (Min)	100 g

Technical data - Construction

Material weighing plate	Stainless steel
Weighing surface (d)	248×198 mm
Dimensions of weighing plate (W×D)	248×198 mm
Material casing	Stainless steel

Technical data - Display

Display height	25 mm
----------------	-------

Technical data - Power supply

Battery	4×15 V AA
Operating time	48 h
Mains adapter/adaptor included	EURO
Input Voltage	230 V AC, 50 Hz

Technical data - Measuring system

Calibration method	Adjusting with external weight
Linearity	± 10 g
Recommended adjusting weight	10 kg (M1)
Repeatability	5 g
Weighing units	Kg
Number of weighing ranges	1
Weighing system	Strain gauge

Technical data - Environmental conditions

Maximum humidity	95 %
Maximum operating temperature	40 °C
Minimum ambient temperature	0 °C
Minimum humidity	25 %

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	420×305×180 mm
------------------------------	----------------

KERN FOB 10K-3LM

Compact stainless steel bench scale with special checkweighing display for even more efficient working, verification optional



Gross weight 4,8 kg

Net weight 3,8 kg

Shipping method Parcel service

Packaging component - by weight - paper 100 g

Packaging component - by weight - cardboard 450 g

Packaging component - by weight - plastics 450 g

Packaging component - by weight - styrofoam 0 g

Packaging component - by weight - wood 0 g

Packaging component - by weight - foam 0 g

Delivery time 1 d

Technical data - Verification

Verification by KERN possible Yes

Verification class III

Optional Services

DAkkS Calibration 963-128

Verification 965-228

Functions

Standard



Option

DAkkS

Factory

M