

## PRODUCT BROCHURE

---

# KERN 308-34

Date: 02.05.2026

ENGLISH



Milligram weights, E1, 10 mg, Wire, stainless steel polished

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 308-34                         |
| Model-Serie           | 308                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML E1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 404576 148 1008                |
| Predecessor Model     | 308-34-2024a                   |
| Customs tariff number | 84239090                       |
| Country of Origin     | CN                             |

### Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 10 mg       |
| Tolerance (OIML) | 0,003 mg    |
| Case aluminium   | 317-009-600 |
| Case plastic     | 347-009-400 |
| Case wood        | 338-090-200 |

### Technical data - Construction

|                    |                          |
|--------------------|--------------------------|
| Dimensions (W×D×H) | 7×7×0,5 mm               |
| Material           | stainless steel polished |

### Technical data - Measuring system

|               |       |
|---------------|-------|
| Nominal value | 10 mg |
|---------------|-------|

|                              |          |
|------------------------------|----------|
| OIML class                   | E1       |
| OIML conform                 | yes      |
| Tolerance (OIML)             | 0,003 mg |
| Density (kg/m <sup>3</sup> ) | 8000     |
| Volume (cm <sup>3</sup> )    | 0,0013   |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 7×7×3 mm       |
| Gross weight                               | 0.05 kg        |
| Net weight                                 | 0.05 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

### Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-254 |
|-------------------|---------|

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|



|        |
|--------|
| Option |
|--------|



## ACCESSORIES

| Model                              | Description                                    |
|------------------------------------|------------------------------------------------|
| <a href="#"><u>338-090-200</u></a> | Wooden Box for Set of Weights KERN 338-090-200 |
| <a href="#"><u>347-009-400</u></a> | Plastic Weight Box KERN 347-009-400            |
| <a href="#"><u>317-009-600</u></a> | Aluminium Protective Box KERN 317-009-600      |
| <a href="#"><u>315-242</u></a>     | Tweezers KERN 315-242                          |
| <a href="#"><u>318-270</u></a>     | Dusting Brush KERN 318-270                     |

## SERVICES

| Model                           | Description                       |
|---------------------------------|-----------------------------------|
| <a href="#"><u>962-254</u></a>  | DAkkS Calibration KERN 962-254    |
| <a href="#"><u>962-254R</u></a> | DAkkS Recalibration KERN 962-254R |

## RELATED ARTICLES DATA

| Model name                     | Nominal value | Tolerance (OIML) | Case aluminium | Case plastic | Case wood   | Weight set case   |
|--------------------------------|---------------|------------------|----------------|--------------|-------------|-------------------|
| <a href="#"><u>308-01</u></a>  | 1 mg          | 0,003 mg         | 317-009-400    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-02</u></a>  | 2 mg          | 0,003 mg         | 317-009-400    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-03</u></a>  | 5 mg          | 0,003 mg         | 317-009-400    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-04</u></a>  | 10 mg         | 0,003 mg         | 317-009-400    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-05</u></a>  | 20 mg         | 0,003 mg         | 317-009-400    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-06</u></a>  | 50 mg         | 0,004 mg         | 317-009-400    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-07</u></a>  | 100 mg        | 0,005 mg         | 317-009-400    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-08</u></a>  | 200 mg        | 0,006 mg         | 317-009-400    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-09</u></a>  | 500 mg        | 0,008 mg         | 317-009-400    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-426</u></a> | 1 mg - 500 mg | -                | -              | -            | -           | in aluminium case |
| <a href="#"><u>308-31</u></a>  | 1 mg          | 0,003 mg         | 317-009-600    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-32</u></a>  | 2 mg          | 0,003 mg         | 317-009-600    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-33</u></a>  | 5 mg          | 0,003 mg         | 317-009-600    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-39</u></a>  | 500 mg        | 0,008 mg         | 317-009-600    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-38</u></a>  | 200 mg        | 0,006 mg         | 317-009-600    | 347-009-400  | 338-090-200 | -                 |
| <a href="#"><u>308-37</u></a>  | 100 mg        | 0,005 mg         | 317-009-600    | 347-009-400  | 338-090-200 | -                 |

| Model name           | Nominal value | Tolerance (OIML) | Case aluminium | Case plastic | Case wood   | Weight set case |
|----------------------|---------------|------------------|----------------|--------------|-------------|-----------------|
| <b><u>308-36</u></b> | 50 mg         | 0,004 mg         | 317-009-600    | 347-009-400  | 338-090-200 | -               |
| <b><u>308-35</u></b> | 20 mg         | 0,003 mg         | 317-009-600    | 347-009-400  | 338-090-200 | -               |
| <b><u>308-34</u></b> | 10 mg         | 0,003 mg         | 317-009-600    | 347-009-400  | 338-090-200 | -               |
| <b><u>308-42</u></b> | 1 mg - 500 mg | -                | -              | 328-222-420  | -           | in plastic case |
| <b><u>308-22</u></b> | 1 mg - 500 mg | -                | -              | -            | -           | in plastic case |