

PRODUCT BROCHURE

KERN YKC-01

Date: 19.09.2025

ENGLISH



Reliable thermal printer for standard applications

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- Affordable thermal printer for printing weights, recipes, quantities, etc. or for GLP recording of weighing data
- Speed: 300 mm/s
- Max. 58 characters per line
- Supported paper widths: 80 mm (roll)
- Operation via mains adapter 230 V/50 Hz (standard EU), as standard
- Data interface RS-232 as standard
- Interface cable not included with delivery

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YKC-01
Model-Serie	YKC
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Main item;Accessories/Services

EAN code (Model)	4045761230378
Predecessor Model	YKC-01-2023a
Customs tariff number	84239010
CE Label	yes

Technical data - Power supply

Mains adapter/adaptor included	EURO
Input Voltage	100 V - 240 V
Mains adapter	Power supply unit

Technical data - Environmental conditions


Maximum humidity	85 %
Maximum operating temperature	50 °C
Minimum ambient temperature	0 °C
Minimum humidity	20 %

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	220×215×150 mm
Gross weight	1.4 kg
Net weight	0.80 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	320 g
Packaging component - by weight - styrofoam	320 g
Delivery time	1 d

FUNCTIONS

Standard



ACCESSORIES

Model	Description
822-903	Power supply KERN 822-903
YKC-A01	Paper Roll KERN YKC-A01
YKC-A10	Paper Rolls KERN YKC-A10

SERVICES

Model	Description
970-015	Warranty extension KERN (+3 years) 970-015