

PRODUCT BROCHURE

KERN 307-13

Date: 09.03.2026

ENGLISH



Individual weight, OIML Class E1, 5 kg, Knob, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class E1, 5 kg, Knob, Stainless steel polished
- Note: For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	307-13
Model-Serie	307
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761002777
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	2,5 mg
Case aluminium	317-130-600
Case plastic	317-130-400
Case wood	317-130-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	5 kg
OIML class	E1
OIML conform	yes
Tolerance (OIML)	2,5 mg
Density (kg/m ³)	8000
Volume (cm ³)	625

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	80×80×140 mm
Gross weight	5 kg
Net weight	5 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	963-243
-------------------	---------

FUNCTIONS

Standard



1 DAY

Option



ACCESSORIES

Model	Description
<u>317-280</u>	Gloves KERN 317-280
<u>317-130-100</u>	Wooden Box for Set of Weights KERN 317-130-100
<u>317-130-400</u>	Plastic Weight Box KERN 317-130-400
<u>317-130-600</u>	Aluminium Protective Box KERN 317-130-600

SERVICES

Model	Description
<u>963-243</u>	DAkkS Calibration KERN 963-243
<u>962-243R</u>	DAkkS Recalibration KERN 962-243R
<u>921-243</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-243