

KERN 333-096

Weight set, F2, 1 mg - 10 kg, Knob, stainless steel fine turned, in aluminium case



Article codes, product groups & Status information

Model article-code	333-096
Model-Serie	333
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main article
EAN code (Model)	4045761224490
Planned start date	2019-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 10 kg
Case aluminium	313-090-600
Weight set case	In aluminium case

Technical data - Construction

Material	Stainless steel fine turned
Is weight set	Yes
Type of case	In aluminium case

Technical data - Measuring system

Nominal value	1 mg - 10 kg
OIML class	F2
OIML conform	Yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	610×260×235 mm
Gross weight	28 kg
Net weight	24 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	4,1 kg
Delivery time	1 d

Optional Services

DAkkS Calibration	962-408
-------------------	---------

Functions

Standard

KERN 333-096

Weight set, F2, 1 mg - 10 kg, Knob, stainless steel fine turned, in aluminium case



Option

