

PRODUCT BROCHURE

KERN 354-076

Date: 10.03.2026

ENGLISH



Weight set, M2, 1 g - 2 kg, Knob, stainless steel fine turned, in aluminium case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	354-076
Model-Serie	354
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main article
EAN code (Model)	404576 1292406
Predecessor Model	354-476
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 2 kg
Case aluminium	314-070-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 g - 2 kg
OIML class	M2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	278×180×138 mm
Gross weight	6 kg
Net weight	6 kg
Shipping method	Parcel service

Packaging component - 400 g
by weight - plastics

Delivery time 1 d

Optional Services

DAkkS Calibration 962-620

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-620	DAkkS Calibration KERN 962-620
962-620R	DAkkS Recalibration KERN 962-620R
920-715	Marking (Last digits of serial number) KERN 920-715

SPARE PARTS

Model	Description
314-070-600	Aluminium Protective Box KERN 314-070-600
315-242	Tweezers KERN 315-242