

PRODUCT BROCHURE

OPTICS OZB-A5213

Date: 30.04.2025

ENGLISH



Individuality, variety and flexible working through our modular construction system ? PREMIUM-Universal stands

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

Stereomicroscope Stand OZB-A5213

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OZB-A5213
Model-Serie	OZB-UP
Assortment	Optics
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761179295
Customs tariff number	90119000

Technical data - Construction

Stand type	universal (ball-bearing double telescopic arm with clamp)
Holder diameter	76 mm

Dimensions (W×D×H)	80×545×515 mm
--------------------	---------------

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	680×330×260 mm
Gross weight	11 kg
Net weight	9 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	1140.53 g
Packaging component - by weight - plastics	22.9164 g

Packaging component - 599.352 g
by weight - styrofoam

Packaging component - 0 g
by weight - foam

Delivery time 1 d

FUNCTIONS

Standard



ACCESSORIES

Model	Description
OZB-A5301	Stereo Microscope Modular System – holders KERN OZB-A5301
OZB-A5306	Stereo Microscope Modular System – holders KERN OZB-A5306

SPARE PARTS

Model	Description
OZB-A5806	Replacement lever OZB-A5806