

KERN 326-05

Individual weight, OIML Class F1, 20 g, Eco Shape, stainless steel polished



Article codes, product groups & Status information

Model article-code	326-05
Model-Serie	326
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761223707
Predecessor Model	326-05-2018a; 326-05K+C; 327-65; 329-05
Successor Model	326-05-2025e
Planned start date	2019-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	20 g
Tolerance (OIML)	0,25 mg
Case aluminium	317-050-600
Case plastic	347-050-400
Case wood	317-050-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	20 g
OIML class	F1

Functions

Standard

OIML conform	Yes
Tolerance (OIML)	0,25 mg
Density (kg/m ³)	7950

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	13×13×22 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	20 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-435
-------------------	---------

KERN 326-05

Individual weight, OIML Class F1, 20 g, Eco Shape, stainless steel polished

