

PRODUCT BROCHURE

KERN 313-054

Date: 09.03.2026

ENGLISH



Weight set, E2, 1 mg - 500 g, Knob, stainless steel polished, in plastic case


CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information



ENGLISH

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

Model article-code	313-054
Model-Serie	313
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main article
EAN code (Model)	4045761081000
Customs tariff number	84239090
Country of Origin	CN

Nominal value	1 mg - 500 g
---------------	--------------

Case plastic	313-052-400
Weight set case	in plastic case

Material	stainless steel polished
Is weight set	yes
Number of weights in set	24
Type of case	in plastic case

Nominal value	1 mg - 500 g
OIML class	E2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	250×210×75 mm
Gross weight	1.6 kg
Net weight	1.6 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-304
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-304</u>	DAkkS Calibration KERN 962-304
<u>962-304R</u>	DAkkS Recalibration KERN 962-304R
<u>921-304</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-304

Model	Description
<u>920-705</u>	Marking (Last digits of serial number) KERN 920-705

SPARE PARTS

Model	Description
<u>313-052-400</u>	Plastic Weight Box KERN 313-052-400