

**KERN
YKA-40**

Date: 04.03.2025

ENGLISH



Image not
available

GB mains adapter, cable length

CONTACT

 Address

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

 Phone and fax

Phone : +49 - [0] 7433 - 9933 - 0
Fax : +49 - [0] 7433 - 9933 - 149

 Internet

Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Power adapter YKA-40 (9 V, 600 mA | 100 - 240 V AC 50/60 Hz 0.5 A | EURO, UK, CH)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YKA-40
Model-Serie	YKA
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Spare part;Accessories/Services
EAN code (Model)	404576 125904 1
Predecessor Model	YKA-40-2020a
Customs tariff number	85044095
CE Label	yes

Packaging component - by weight - cardboard	0 g
---	-----

Packaging component - by weight - plastics	0 g
--	-----

Packaging component - by weight - glass	0 g
---	-----

Packaging component - by weight - wood	0 g
--	-----

Delivery time	1 d
---------------	-----

ENGLISH

Technical data - Power supply

Mains adapter/adapter included	CH;EURO;UK
Input Voltage	100-240V AC 50/60Hz 0.5A

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	93×67×70 mm
Gross weight	0.20 kg
Net weight	0.15 kg
Shipping method	Parcel service

ACCESSORIES

Model	Description
YKA-40-US	Plug SAUTER YKA-40-US