

OZB-A5301



Holders

Features

- There are two microscope head holders available for these flexible, modular systems. These brackets are suitable for all stereo microscope stands and universal stands (excluding spring loaded arm), to make focusing possible
- The first variant available is a holder with adjustable handwheel as well as adjustment of the torque for your configuration
- For professional applications you have the choice of a mount with coarse and fine focusing knob for the very best focusing operation
- Diameter of the connector for the microscope head: 76 mm
- Diameter of the connector for the stand:25 mm

| Model | Description | |
|-----------|--|--|
| KERN | | |
| OZB-A5301 | Holder with adjustable torque of the hand wheel. Suitable for all universal stands (except of spring loaded arm) and for all basic stands as possible accessories. | |
| OZB-A5306 | Holder with coaxial coarse and fine focusing knob and adjustable torque of the hand wheel. Suitable for all universal stands (except of spring loaded arm) and for all basic stands as possible accessories. | |

Stereomicroscope Modular System - Dust Covers KERN OBB-C



Dust covers

Features

- When working with microscopes, we offer dust covers to give greater ease of use. By using these, you can easily avoid the time-consuming cleaning work which is necessary with routine use of your microscope
- Depending on the size of your microscope set or your microscope configuration you can select between three different models
- Please find detailed information in the following model outfit list

| Model | Description | Suitable for |
|-----------|--------------------|--|
| KERN | | |
| OBB-A1387 | Size 1: 485×450 mm | Stereo microscope heads |
| OBB-A1388 | Size 2: 600×650 mm | Stereo microscope heads in combination with basic stands |
| OBB-A1389 | Size 3: 700×900 mm | Stereo microscope sets, stereo microscope heads in combination with universal stands |