

## PRODUCT BROCHURE

---

# KERN 348-08

Date: 10.03.2026

ENGLISH



Milligram weights, M1, 200 mg, Platelet, stainless steel polished

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

- Milligram weights, M1, 200 mg, Platelet, stainless steel polished (OIML)
- For the verification you need a box. Please add a box for single weights.

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 348-08                         |
| Model-Serie           | 348                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML M1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761240209                  |
| Predecessor Model     | 348-08-2022a                   |
| Customs tariff number | 84239090                       |
| Country of Origin     | CN                             |

### Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 200 mg      |
| Tolerance (OIML) | 0,6 mg      |
| Case aluminium   | 317-009-600 |
| Case plastic     | 347-009-400 |
| Case wood        | 338-090-200 |

### Technical data - Construction

|                    |          |
|--------------------|----------|
| Dimensions (W×D×H) | 9×9×4 mm |
|--------------------|----------|

Material stainless steel polished

### Technical data - Measuring system

|                              |        |
|------------------------------|--------|
| Nominal value                | 200 mg |
| OIML class                   | M1     |
| OIML conform                 | yes    |
| Tolerance (OIML)             | 0,6 mg |
| Density (kg/m <sup>3</sup> ) | 7960   |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 9×9×4 mm       |
| Gross weight                               | 0.05 kg        |
| Net weight                                 | 0.05 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

### Technical data - Verification

|                       |     |
|-----------------------|-----|
| Verification possible | yes |
|-----------------------|-----|

### Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-658 |
|-------------------|---------|

## FUNCTIONS

Standard



Option



ENGLISH

## ACCESSORIES

| Model                              | Description                                    |
|------------------------------------|------------------------------------------------|
| <a href="#"><u>338-090-200</u></a> | Wooden Box for Set of Weights KERN 338-090-200 |
| <a href="#"><u>347-009-400</u></a> | Plastic Weight Box KERN 347-009-400            |
| <a href="#"><u>317-009-600</u></a> | Aluminium Protective Box KERN 317-009-600      |
| <a href="#"><u>315-242</u></a>     | Tweezers KERN 315-242                          |
| <a href="#"><u>318-270</u></a>     | Dusting Brush KERN 318-270                     |

## SERVICES

| Model                           | Description                       |
|---------------------------------|-----------------------------------|
| <a href="#"><u>962-658</u></a>  | DAkks Calibration KERN 962-658    |
| <a href="#"><u>962-658R</u></a> | DAkks Recalibration KERN 962-658R |