

PRODUCT BROCHURE

SAUTER AFM 15

Date: 30.04.2025

ENGLISH



Inner threaded sleeve, steel; M12 threaded hole - M6 threaded hole

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

Adapter (Threaded hole M12, Threaded hole M6) AFM 15

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	AFM 15
Model-Serie	AFM
Assortment	Sauter
Product Type	Force measurement
Product group	Test stands (motorised)
Product usage type	Accessories/DL
EAN code (Model)	404576 1244955
Customs tariff number	90278990

Technical data - Construction

Mounting threads	Threaded hole M12
------------------	-------------------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	30×15×15 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
UN number	No dangerous goods
Packaging component - by weight - plastics	0 g
Delivery time	1 d

FUNCTIONS

Standard

