

EU Declaration of Conformity | EU-Konformitätserklärung

EN Wir erklären hiermit unter alleiniger Verantwortung, dass das Produkt, auf das sich diese Erklärung bezieht, mit den nachstehenden Richtlinien übereinstimmt. Für die Waage liegt eine EG-Bauartzulassung vor. Das Jahr der EG-Eichung ist neben der CE Konformitätskennzeichnung aufgeführt. Solche Waagen sind ab Werk EG-geeicht und tragen die Kennzeichnung „M“ auf dem Gerät selbst und auf der Verpackung. Sie dürfen im eichpflichtigen Verkehr verwendet werden. Die Fallbeschleunigung gibt bei EG-geeichten Waagen an, für welchen Aufstellungsort die Waage EG-geeicht ist.

DE We hereby declare and assume sole responsibility for the declaration that the product complies with the standards hereinafter. The scale was issued with an EC Type Approval Certificate. The year of EC verification can be found next to the CE mark of conformity. Such scales are EC verified ex factory and bear the mark “M”, on the device itself as well as on the packaging. They can be used for weighing with obligation to verify. For EC verified scales, gravity acceleration indicates for which place of installation the scale is EC verified.

Model | Modell |

Type | Typ |


KDP 10K-3
KDP 10K-4
KDP 3000-2
KDP 300-3

TKDP 10K-3-A
TKDP 10K-4-A
TKDP 3000-2-A
TKDP 300-3-A

CE mark applied CE Kennzeichnung	EU directive EU-Richtlinie	Standards Normen
	2011/65/EU (RoHS)	EN 63000:2018
	2014/30/EU (EMC)	EN 61326-1:2013 EN 61326-2-2:2013
	2014/35/EU (LVD)	EN 61010-1:2010 EN 62368-1:2014+A11:2017
	2014/53/EU (RED)	EN 62311:2008 ETSI EN 301 489- 1 V2.2.0 (draft, 2017-03) ETSI EN 301 489- 17 V3.2.0 (draft, 2017-03)

Date | Datum | : 10.11.2021

Place of issue: 72336 Balingen,
Ort der Ausstellung: Germany



Albert Sauter
KERN & Sohn GmbH

Signature: Managing director
Signatur: Geschäftsführer