


# KERN 346-042-008

Weight container for rectangular weights, calibrated to OIML class M1



## Article codes, product groups & Status information

|                       |                      |
|-----------------------|----------------------|
| Model article-code    | 346-042-008          |
| Model-Serie           | 346                  |
| Assortment            | KERN                 |
| Product Type          | Test weights         |
| Product group         | OIML M1 Test weights |
| Product usage type    | Main article         |
| EAN code (Model)      | 4045761225060        |
| Planned start date    | 2019-1HJ             |
| Customs tariff number | 84239090             |

|                                            |         |
|--------------------------------------------|---------|
| Net weight                                 | 40 kg   |
| Shipping method                            | Freight |
| Packaging component - by weight - plastics | 0 g     |
| Delivery time                              | 2 d     |

## Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-646 |
|-------------------|---------|

## Technical data - Primary

|                  |       |
|------------------|-------|
| Nominal value    | 40 kg |
| Tolerance (OIML) | 1,5 g |

## Technical data - Construction

|                    |                 |
|--------------------|-----------------|
| Dimensions (W×D×H) | 710×300×470 mm  |
| Material           | Stainless steel |

## Technical data - Measuring system

|                  |        |
|------------------|--------|
| Nominal value    | 40 kg  |
| OIML class       | M1     |
| OIML conform     | Yes    |
| Tolerance (OIML) | 1,5 mg |

## Technical data - Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Dimensions packaging (W×D×H) | 710×300×470 mm |
| Gross weight                 | 40 kg          |

## Functions

# KERN 346-042-008

Weight container for rectangular weights, calibrated to OIML class M1



Standard

