

## PRODUCT BROCHURE

---

# KERN 315-207-204

Date: 19.03.2025

ENGLISH



Safety Set 10 g (F1); 100 g (F1)

---

## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)

## PRODUCT SHOWCASE



## DESCRIPTION

- These “KERN Safety Sets” each contain one test weight for testing the sensitivity, i.e. the correct adjustment of your balance as well as a small test weight to test the lower end of the weighing range, the so-called minimum sample weight
- As an option, the “KERN Safety Set” has space for another test weight, for testing your balance at a weight which is relevant for you
- Useful accessories which have been selected to suit that particular “KERN Safety Set”, such as, for example, special gloves, tweezers, weight grips, brushes, etc. will assist you in handling your test weights properly
- Stored in the practical protective case next to your balance, you can check and ensure the high precision of your balance at any time

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	315-207-204
Model-Serie	315
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services

EAN code (Model) 4045761245471

Customs tariff number 84239090

### Technical data - Primary

Nominal value 10 g & 100 g

Case aluminium zet315-000-000

Weight set case in aluminium case

### Technical data - Construction

Material stainless steel

Is weight set	yes
Number of weights in set	2
Type of case	in aluminium case

## Technical data - Measuring system

Nominal value	10 g & 100 g
OIML class	F1
OIML conform	yes

## Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	150×115×70 mm
Gross weight	1.00 kg
Net weight	0.50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	500 g
Delivery time	1 d

## Technical data - Verification

Verification possible	yes
-----------------------	-----

## Optional Services

DAkkS Calibration	962-400-400
-------------------	-------------

## FUNCTIONS

Standard
----------



## SERVICES

Model	Description
962-400-400	DAkkS Calibration KERN 962-400-400

## SPARE PARTS

Model	Description
315-242	Tweezers KERN 315-242

## RELATED ARTICLES DATA

Model name	Nominal value	Case aluminium	Weight set case
315-208-205	20 g & 200 g	zet315-000-000	in aluminium case
315-207-204	10 g & 100 g	zet315-000-000	in aluminium case