

KERN 347-057

Individual weight, OIML Class M1, 20 g, Slot weight (Type B), stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-057
Model-Serie	347-057
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366657
Planned start date	2024-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	20 g
Tolerance (OIML)	2,5 mg

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	20 g
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	2,5 mg
Density (kg/m ³)	8.000
Volume (cm ³)	2,381

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	38×38×4 mm
------------------------------	------------

Gross weight	0,05 kg
--------------	---------

Net weight	0,05 kg
------------	---------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - plastics	0 g
--	-----

Delivery time	1 d
---------------	-----

Optional Services

DAkkS Calibration	962-635
-------------------	---------

Functions

KERN 347-057

Individual weight, OIML Class M1, 20 g, Slot weight (Type B), stainless steel fine turned



Standard



Option

