

KERN 367-11

Individual weight, OIML Class M3, 1 kg, Knob, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	367-11
Model-Serie	367-11
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761257290
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 kg
Tolerance (OIML)	500 mg
Case aluminium	317-110-600
Case plastic	347-110-400
Case wood	337-110-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	1 kg
OIML class	M3
OIML conform	Yes
Tolerance (OIML)	500 mg

Functions

Standard

Density (kg/m ³)	7.950
Volume (cm ³)	119,0476

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	48×48×80 mm
Gross weight	1,00 kg
Net weight	1,00 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	No
-----------------------	----

Optional Services

DAkkS Calibration	962-641
-------------------	---------

KERN 367-11

Individual weight, OIML Class M3, 1 kg, Knob, stainless steel fine turned



Option

