

# KERN 328-35

Milligram weights, F1, 20 mg, Wire, stainless steel polished



## Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 328-35                         |
| Model-Serie           | 328                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML F1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761478626                  |
| Predecessor Model     | 328-35-2.024a                  |
| Planned start date    | 2025-2HJ                       |
| Customs tariff number | 84239090                       |
| Country of Origin     | CN                             |

## Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 20 mg       |
| Tolerance (OIML) | 0,03 mg     |
| Case aluminium   | 317-009-600 |
| Case plastic     | 347-009-400 |
| Case wood        | 338-090-200 |

## Technical data - Construction

|                    |                          |
|--------------------|--------------------------|
| Dimensions (W×D×H) | 7×7×3 mm                 |
| Material           | Stainless steel polished |

## Technical data - Measuring system

|               |       |
|---------------|-------|
| Nominal value | 20 mg |
| OIML class    | F1    |

## Functions

Standard

|                              |         |
|------------------------------|---------|
| OIML conform                 | Yes     |
| Tolerance (OIML)             | 0,03 mg |
| Density (kg/m <sup>3</sup> ) | 8.000   |
| Volume (cm <sup>3</sup> )    | 0,0025  |

## Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 7×7×3 mm       |
| Gross weight                               | 0,05 kg        |
| Net weight                                 | 0,05 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

## Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-455 |
|-------------------|---------|

# KERN 328-35

Milligram weights, F1, 20 mg, Wire, stainless steel polished



Option

