

PRODUCT BROCHURE

KERN 315-208-204

Date: 19.03.2025

ENGLISH



Safety Set 10 g (F1); 200 g (F1)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



DESCRIPTION

- These “KERN Safety Sets” each contain one test weight for testing the sensitivity, i.e. the correct adjustment of your balance as well as a small test weight to test the lower end of the weighing range, the so-called minimum sample weight
- As an option, the “KERN Safety Set” has space for another test weight, for testing your balance at a weight which is relevant for you
- Useful accessories which have been selected to suit that particular “KERN Safety Set”, such as, for example, special gloves, tweezers, weight grips, brushes, etc. will assist you in handling your test weights properly
- Stored in the practical protective case next to your balance, you can check and ensure the high precision of your balance at any time

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	315-208-204
Model-Serie	315
Assortment	KERN
Product Type	Test weights
Product group	Safety sets
Product usage type	Main item;Accessories/Services

EAN code (Model)	4045761245693
------------------	---------------

Customs tariff number	84239090
-----------------------	----------

Technical data - Measuring system

Nominal weights, included	10 g (F1); 200 g (F1)
---------------------------	-----------------------

Nominal value	10 g & 200 g
---------------	--------------

OIML class	F1
------------	----

OIML conform	yes
--------------	-----

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	150×115×70 mm
Gross weight	1.00 kg
Net weight	0.50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	500 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-400-400
-------------------	-------------

FUNCTIONS

Standard



SERVICES

Model	Description
962-400-400	DAkkS Calibration KERN 962-400-400

SPARE PARTS

Model	Description
315-242	Tweezers KERN 315-242