

KERN OZB-A4572

Professional illumination units guarantee outstanding, uniform and strong illumination



Article codes, product groups & Status information

Model article-code	OZB-A4572
Model-Serie	OZB-IR
Assortment	Optics
Product Type	Microscopes
Product group	External Illumination
Product usage type	Main article
EAN code (Model)	4045761218437
Predecessor Model	OZB-A4572-2018a
Planned start date	2019-1HJ
Customs tariff number	94054239

Technical data - Primary

Illumination type	LED
Illumination intensity	4W

Technical data - Construction

Color tempertaure	6500-7000 K
Dimmable	Yes
By segments	Yes
Ring diameter	6 cm

Technical data - Power supply

Protection Class	Class III - Low voltage
------------------	-------------------------

Technical data - Environmental conditions

Storage Temp Min	-5 °C
------------------	-------

Functions

Standard

Stotage Temp Max	40 °C
------------------	-------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	235×200×50 mm
Gross weight	0.50 kg
Net weight	0.45 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	72.0324 g
Packaging component - by weight - plastics	2.3112 g
Packaging component - by weight - styrofoam	21.9564 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Verification

WEEE	Yes
------	-----

KERN OZB-A4572

Professional illumination units guarantee outstanding, uniform and strong illumination

