

KERN 328-33

Milligram weights, F1, 5 mg, Wire, stainless steel polished



Article codes, product groups & Status information

Model article-code	328-33
Model-Serie	328
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761478602
Predecessor Model	328-33-2.024a
Planned start date	2025-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	5 mg
Tolerance (OIML)	0,02 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	12×12×3 mm
Material	Stainless steel polished

Technical data - Measuring system

Nominal value	5 mg
OIML class	F1

Functions

Standard

OIML conform Yes

Tolerance (OIML) 0,02 mg

Density (kg/m³) 8.000

Volume (cm³) 0,0006

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 12×12×3 mm

Gross weight 0,05 kg

Net weight 0,05 kg

Shipping method Parcel service

Packaging component - by weight - plastics 0 g

Delivery time 1 d

Optional Services

DAkkS Calibration 962-453

KERN 328-33

Milligram weights, F1, 5 mg, Wire, stainless steel polished



Option

