


# KERN 348-06

Milligram weights, M1, 50 mg, Platelet, stainless steel polished



## Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 348-06                         |
| Model-Serie           | 348                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML M1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761240186                  |
| Predecessor Model     | 348-6-2.022a                   |
| Planned start date    | 2020-1HJ                       |
| Customs tariff number | 84239090                       |
| Country of Origin     | CN                             |

## Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 50 mg       |
| Tolerance (OIML) | 0,4 mg      |
| Case aluminium   | 317-009-600 |
| Case plastic     | 347-009-400 |
| Case wood        | 338-090-200 |

## Technical data - Construction

|                    |                          |
|--------------------|--------------------------|
| Dimensions (W×D×H) | 7,5×6×3 mm               |
| Material           | Stainless steel polished |

## Technical data - Measuring system

|               |       |
|---------------|-------|
| Nominal value | 50 mg |
| OIML class    | M1    |

## Functions

Standard

|                              |        |
|------------------------------|--------|
| OIML conform                 | Yes    |
| Tolerance (OIML)             | 0,4 mg |
| Density (kg/m <sup>3</sup> ) | 7.960  |

## Technical data - Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Dimensions packaging (W×D×H) | 7,5×6×3 mm     |
| Gross weight                 | 0,05 kg        |
| Net weight                   | 0,05 kg        |
| Shipping method              | Parcel service |
| Delivery time                | 1 d            |

## Optional Services

|                   |         |
|-------------------|---------|
| DAkKS Calibration | 962-656 |
|-------------------|---------|

# KERN 348-06

Milligram weights, M1, 50 mg, Platelet, stainless steel polished



Option

