

PRODUCT BROCHURE

KERN 323-09

Date: 09.03.2026

ENGLISH



Weight set, F1, 1 mg - 10 kg, Knob, stainless steel polished, in wooden case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	323-09
Model-Serie	323
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1004320
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 10 kg
---------------	--------------

Case wood	313-090-100
Weight set case	in wooden case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Number of weights in set	29
Type of case	in wooden case

Technical data - Measuring system

Nominal value	1 mg - 10 kg
OIML class	F1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	540×200×240 mm
Gross weight	34 kg
Net weight	32 kg
Shipping method	Freight
Packaging component - by weight - plastics	1890 g
Delivery time	2 d

Optional Services

DAkkS Calibration	962-408
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-408</u>	DAkkS Calibration KERN 962-408
<u>962-408R</u>	DAkkS Recalibration KERN 962-408R
<u>921-408</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-408
<u>920-709</u>	Marking (Last digits of serial number) KERN 920-709

SPARE PARTS

Model	Description
<u>313-090-100</u>	Box of Weights KERN 313-090-100
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270