

## PRODUCT BROCHURE

---

# KERN 312-024

Date: 09.03.2026

ENGLISH



Weight set, E2, 1 g - 50 g, Compact form, stainless steel polished, in plastic case

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                      |
|-----------------------|----------------------|
| Model article-code    | 312-024              |
| Model-Serie           | 312                  |
| Assortment            | KERN                 |
| Product Type          | Test weights         |
| Product group         | OIML E2 Test weights |
| Product usage type    | Main article         |
| EAN code (Model)      | 4045761116559        |
| Customs tariff number | 84239090             |
| Country of Origin     | DE                   |

### Technical data - Primary

|                 |                 |
|-----------------|-----------------|
| Nominal value   | 1 g - 50 g      |
| Case plastic    | 313-052-400     |
| Weight set case | in plastic case |

### Technical data - Construction

|                          |                          |
|--------------------------|--------------------------|
| Material                 | stainless steel polished |
| Is weight set            | yes                      |
| Number of weights in set | 8                        |
| Type of case             | in plastic case          |

### Technical data - Measuring system

|               |            |
|---------------|------------|
| Nominal value | 1 g - 50 g |
| OIML class    | E2         |
| OIML conform  | yes        |

### Technical data - Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Dimensions packaging (W×D×H) | 250×210×75 mm  |
| Gross weight                 | 0.50 kg        |
| Net weight                   | 0.50 kg        |
| Shipping method              | Parcel service |

Packaging component - 0 g  
by weight - plastics

Delivery time 1 d

## Technical data - Verification

Verification possible no

## Optional Services

DAkkS Calibration 962-315

## FUNCTIONS

Standard



Option



## SERVICES

| Model                    | Description                                          |
|--------------------------|------------------------------------------------------|
| <a href="#">962-315</a>  | DAkkS Calibration KERN 962-315                       |
| <a href="#">962-315R</a> | DAkkS Recalibration KERN 962-315R                    |
| <a href="#">921-315</a>  | Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-315 |
| <a href="#">920-710</a>  | Marking (Last digits of serial number) KERN 920-710  |

## SPARE PARTS

| Model                       | Description                         |
|-----------------------------|-------------------------------------|
| <a href="#">313-052-400</a> | Plastic Weight Box KERN 313-052-400 |
| <a href="#">315-242</a>     | Tweezers KERN 315-242               |