

KERN SEQ-01

Equipment Qualification made easy – with the KERN Equipment Qualification Software (EQS)



Equipment qualification **KERN**
CALIBRATION

Logbook for qualification of balance and weighing systems

General details of equipment

Device Type:	Single-range balance: ☒	Multi-range balance: ☐	Multi-interval balance: ☐
Verified Weighing Device: ☐			
	Range 1	Range 2	Range 3
Weighing range (Max):	335 g	N/A	N/A
Scale interval (d):	0.0001 g	N/A	N/A
Verification value (e):	0.001 g	N/A	N/A
Min:	0.01 g	N/A	N/A
Connected accessories:	Printer YNG-01		
Operating instructions managed by:	Process owner		
Date of installation:	29-08-2023		
Date of putting out of service:	N/A		
Service contract number:	N/A		
Service case telephone number:	N/A		
Equipment qualification executed by:	Qualification officer		
Time Günther			
(Name of qualification officer)	(signature)		
Equipment qualification reviewed by:	John Doe		
John Doe			
(Name of measuring and test equipment officer)	(signature)		

Article codes, product groups & Status information

Model article-code	SEQ-01
Model-Serie	SEQ
Assortment	Scales
Product Type	Software
Product group	Software (balances), other
Product usage type	Main article
EAN code (Model)	4045761361720
Planned start date	2024-1HJ
Customs tariff number	85234910
Country of Origin	DE

Technical data - Packing & Shipping

Packaging component - by weight - paper	100 g
Packaging component - by weight - cardboard	0 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - wood	0 g
Packaging component - by weight - foam	0 g
Delivery time	1 d