

PRODUCT BROCHURE

KERN 308-36

Date: 02.05.2026

ENGLISH



Milligram weights, E1, 50 mg, Wire, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	308-36
Model-Serie	308
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 148 1015
Predecessor Model	308-36-2024a
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	50 mg
Tolerance (OIML)	0,004 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	12×12×0,5 mm
Material	stainless steel polished

Technical data - Measuring system

Nominal value	50 mg
---------------	-------

OIML class	E1
OIML conform	yes
Tolerance (OIML)	0,004 mg
Density (kg/m ³)	8000
Volume (cm ³)	0,0063

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	12×12×3 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-256
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>338-090-200</u>	Wooden Box for Set of Weights KERN 338-090-200
<u>347-009-400</u>	Plastic Weight Box KERN 347-009-400
<u>317-009-600</u>	Aluminium Protective Box KERN 317-009-600
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270

SERVICES

Model	Description
<u>962-256</u>	DAkkS Calibration KERN 962-256
<u>962-256R</u>	DAkkS Recalibration KERN 962-256R

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood	Weight set case
<u>308-01</u>	1 mg	0,003 mg	317-009-400	347-009-400	338-090-200	-
<u>308-02</u>	2 mg	0,003 mg	317-009-400	347-009-400	338-090-200	-
<u>308-03</u>	5 mg	0,003 mg	317-009-400	347-009-400	338-090-200	-
<u>308-04</u>	10 mg	0,003 mg	317-009-400	347-009-400	338-090-200	-
<u>308-05</u>	20 mg	0,003 mg	317-009-400	347-009-400	338-090-200	-
<u>308-06</u>	50 mg	0,004 mg	317-009-400	347-009-400	338-090-200	-
<u>308-07</u>	100 mg	0,005 mg	317-009-400	347-009-400	338-090-200	-
<u>308-08</u>	200 mg	0,006 mg	317-009-400	347-009-400	338-090-200	-
<u>308-09</u>	500 mg	0,008 mg	317-009-400	347-009-400	338-090-200	-
<u>308-426</u>	1 mg - 500 mg	-	-	-	-	in aluminium case
<u>308-31</u>	1 mg	0,003 mg	317-009-600	347-009-400	338-090-200	-
<u>308-32</u>	2 mg	0,003 mg	317-009-600	347-009-400	338-090-200	-
<u>308-33</u>	5 mg	0,003 mg	317-009-600	347-009-400	338-090-200	-
<u>308-39</u>	500 mg	0,008 mg	317-009-600	347-009-400	338-090-200	-
<u>308-38</u>	200 mg	0,006 mg	317-009-600	347-009-400	338-090-200	-
<u>308-37</u>	100 mg	0,005 mg	317-009-600	347-009-400	338-090-200	-

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood	Weight set case
<u>308-36</u>	50 mg	0,004 mg	317-009-600	347-009-400	338-090-200	-
<u>308-35</u>	20 mg	0,003 mg	317-009-600	347-009-400	338-090-200	-
<u>308-34</u>	10 mg	0,003 mg	317-009-600	347-009-400	338-090-200	-
<u>308-42</u>	1 mg - 500 mg	-	-	328-222-420	-	in plastic case
<u>308-22</u>	1 mg - 500 mg	-	-	-	-	in plastic case