

PRODUCT BROCHURE

KERN 326-036

Date: 19.09.2025

ENGLISH



Weight set, F1, 1 g - 100 g, Eco Shape, stainless steel polished, in aluminium case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	326-036
Model-Serie	326
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1224 131
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 g - 100 g
Case aluminium	314-030-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Is weight set yes

Type of case in aluminium case

Technical data - Measuring system

Nominal value	1 g - 100 g
OIML class	F1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	165×125×75 mm
Gross weight	0.90 kg
Net weight	0.50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	440 g
Delivery time	1 d

Optional Services

DAkkS Calibration 962-416

FUNCTIONS

Standard



SERVICES

Model	Description
962-416	DAkkS Calibration KERN 962-416
962-416R	DAkkS Recalibration KERN 962-416R
921-416	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-416
920-711	Marking (Last digits of serial number) KERN 920-711

SPARE PARTS

Model	Description
314-030-600	Aluminium Protective Box KERN 314-030-600
315-242	Tweezers KERN 315-242

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood	Weight set case
326-36	5 kg	25 mg	346-060-600	not available	not available	-
326-37	10 kg	50 mg	346-070-600	not available	not available	-
326-38	20 kg	100 mg	346-080-600	not available	not available	-
326-39	50 kg	250 mg	346-090-600	not available	not available	-
326-022	1 g - 50 g	-	-	-	313-020-100	in wooden case
326-024	1 g - 50 g	-	-	313-052-400	-	in plastic case
326-026	1 g - 50 g	-	314-020-600	-	-	in aluminium case
326-032	1 g - 100 g	-	-	-	313-030-100	in wooden case
326-034	1 g - 100 g	-	-	313-052-400	-	in plastic case

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood	Weight set case
326-036	1 g - 100 g	-	314-030-600	-	-	in aluminium case
326-042	1 g - 200 g	-	-	-	313-040-100	in wooden case
326-044	1 g - 200 g	-	-	313-052-400	-	in plastic case
326-046	1 g - 200 g	-	314-040-600	-	-	in aluminium case
326-052	1 g - 500 g	-	-	-	313-050-100	in wooden case
326-054	1 g - 500 g	-	-	313-052-400	-	in plastic case
326-056	1 g - 500 g	-	314-050-600	-	-	in aluminium case
326-062	1 g - 1 kg	-	-	-	313-060-100	in wooden case
326-064	1 g - 1 kg	-	-	313-082-400	-	in plastic case
326-066	1 g - 1 kg	-	314-060-600	-	-	in aluminium case
326-072	1 g - 2 kg	-	-	-	313-070-100	in wooden case
326-074	1 g - 2 kg	-	-	313-082-400	-	in plastic case
326-076	1 g - 2 kg	-	314-070-600	-	-	in aluminium case
326-082	1 g - 5 kg	-	-	-	313-080-100	in wooden case
326-084	1 g - 5 kg	-	-	313-082-400	-	in plastic case
326-086	1 g - 5 kg	-	314-080-600	-	-	in aluminium case
326-092	1 g - 10 kg	-	-	-	313-090-100	in wooden case
326-096	1 g - 10 kg	-	314-090-600	-	-	in aluminium case
326-01	1 g	0,1 mg	317-010-600	347-030-400	317-010-100	-
326-02	2 g	0,12 mg	317-020-600	347-030-400	317-020-100	-
326-03	5 g	0,16 mg	317-030-600	347-030-400	317-030-100	-
326-04	10 g	0,2 mg	317-040-600	347-050-400	317-040-100	-
326-05	20 g	0,25 mg	317-050-600	347-050-400	317-050-100	-
326-06	50 g	0,3 mg	317-060-600	347-070-400	317-060-100	-
326-07	100 g	0,5 mg	317-070-600	347-070-400	317-070-100	-

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood	Weight set case
326-08	200 g	1 mg	317-080-600	347-080-400	317-080-100	-
326-09	500 g	2,5 mg	317-090-600	347-090-400	317-090-100	-
326-11	1 kg	5 mg	317-110-600	347-110-400	317-110-100	-
326-12	2 kg	10 mg	317-120-600	347-120-400	317-120-100	-
326-13	5 kg	25 mg	317-130-600	347-130-400	317-130-100	-
326-14	10 kg	50 mg	317-140-600	347-140-400	317-140-100	-